

USB Video-Rhino-Laryngoscope

RSX USB / RSX-P Paediatrics



Video-Rhino-Laryngoscope

| | | |
|-------------------|---|--------------|
| Model: | RSX USB | RSX-P |
| Diameters: | 3.9mm | 2.9mm |
| Working length: | 310mm | |
| FOV: | 90° | |
| Tip Articulation: | 130°/130° | |
| Weight: | approx. 330g incl. USB-Connection cable | |
| Interface: | USB | |
| Lightsource: | Integrated high power LED lightsource | |

Completely immersible machine cleaning/ disinfection



Medical Tablet

| | |
|-------------------|--|
| Model: | 10,1" Medical Tablet PC, fanless design |
| Capacity: | 256GB SSD |
| Operating system: | Windows 10 |
| Power supply: | Exchangeable battery |
| Resolution: | 1280x800 Pixel |
| Operation: | Multi-Touch-Panel reinforced glass front |
| Interfaces: | W-LAN 802.11 a/b/g/n, Bluetooth, USB |
| Speakers: | integrated |
| Microphone: | integrated |
| Dimensions: | 292x196x20mm |
| Weight: | 1.1Kg |

USB Video-Rhino-Laryngoscope




from 2.9mm ∅



RSX USB / RSX-P Paediatrics

Plug & Play: All you need

Efficient Diagnostics. Made in Germany




- RSX USB 3.9mm diameter
- RSX-P 2.9mm diameter (Paediatrics)
- ultra-compact with integrated LED
- USB connection




Configurations




Portable system



directly connection to Medical PC



directly connection to Medical Tablet



Medical Tablet